



AFRI MEDICAL COMPANY

THE RIGHT WAY...
IS THE ONLY WAY...



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

SURGICAL DRESSING (OpiFix)



OpiFix PU Dressing

OpiFix Creates new standards in postoperative care.

OpiFix Offers patient's both comfort and protection.

OpiFix Combines protection, conformability, and transparency of polyurethane film with the cushioning and absorbency of a non-adherent pad.

OpiFix Provides a unique barrier system for an easy application.

OpiFix Is a high waterproof dressing. OpiFix reduce risk of maceration

Indications:

- For surgical wounds.
- After injuries (cuts, grazes, etc....)
- For exudating wounds

Availability

Code	Description	Size	Pack
39720	OpiFix	6 cm x 10 cm	100 Dressing
39721	OpiFix	10 cm x 10 cm	40 Dressing
39722	OpiFix	15 cm x 10 cm	40 Dressing
39723	OpiFix	20 cm x 10 cm	40 Dressing
39724	OpiFix	25 cm x 10 cm	40 Dressing
39725	OpiFix	30 cm x 10 cm	40 Dressing
39726	OpiFix	35 cm x 10 cm	40 Dressing
39727	OpiFix	40 cm x 10 cm	40 Dressing
39728	OpiFix (without Pad)	6 cm x 10 cm	400 Dressing
39730	OpiFix (without Pad)	15 cm x 10 cm	100 Dressing



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

SURGICAL DRESSING (ProFix)



ProFix Dressing

ProFix Is a conformable, soft, postoperative, non-woven dressing that assures efficient wound healing process in a very short time.

ProFix With a non-adherent pad absorbs blood and exudates maintaining a good healing environment.

ProFix Provides safe and secured dressing fixation over the wound.

ProFix Is ideally suited for postoperative use, laceration, grazes and sutured wounds.

ProFix Has an easy - peel pack.

Indications:

- For surgical wounds.
- After injuries (cuts, grazes, etc....)
- For exudating wounds

Availability

Code	Description	Size	Pack
34710	ProFix	6 cm x 10 cm	100 Dressing
34711	ProFix	10 cm x 10 cm	40 Dressing
34712	ProFix	15 cm x 10 cm	40 Dressing
34713	ProFix	20 cm x 10 cm	40 Dressing
34714	ProFix	25 cm x 10 cm	40 Dressing
34715	ProFix	30 cm x 10 cm	40 Dressing
34716	ProFix	35 cm x 10 cm	40 Dressing
34717	ProFix	40 cm x 10 cm	40 Dressing
34720	ProFix I.V. Cannula Fixation	6 cm x 10 cm	150 Dressing
34721	ProFix Eye Patch	10 cm x 10 cm	100 Dressing



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

SURGICAL DRESSING (DraPex)



DraPex Incise Drapes

DraPex Reduces the risk of intra- operative wound infection.

DraPex Is a transparent, thin glaze- free film drape, waterproof and a complete bacterial barrier

DraPex Prevents access of environmental bacteria to the wound, because of the high vapor permeability of the film the drapes stays closely adhered to the skin and wound margin throughout the procedures, locking bacteria into the skin ducts and structures.

DraPex Stretches and regains to allow skin retraction without trauma to the skin surface, they mould to the body contour and permit free manipulation of the joints during orthopedic surgery.

Application of the drape.

- Clean the skin in the normal manner and allow the antibacterial skin prep, to dry slightly 10-15 seconds apply the linen or disposable drape.
- Select a appropriate size of Drapex and remove from packet, holding both green handle edge taught across the patient a third scrapped member of the team removes the packing paper.
- This is best achieved by pulling diagonally from one corner.
- Lower one green handle to the surface of the draped patient and gently press the adhesive into place smoothing the film across the exposed skin lower the film gently into position.
- Before starting to cut, use a blunt instrument to press the film to the skin along the proposed incision line.
- At the end of the surgery, there are two options:
 - 1- Remove the Drapex prior to placing the final skin closures.
 - 2- Close the incision through the film, in this option the film is wiped dry and cut away from linen drape, a postoperative dressing is placed over the suture line treating the film as if it were the skin.

Drapex incise drapes are sterilized by ethylene oxide and are available in the following sizes.

Availability

Code	Description	Size	Pack
37704	DraPex	10 cm x 15 cm	20 Dressing
37705	DraPex	15 cm x 28 cm	20 Dressing
37707	DraPex	40 cm x 42 cm	10 Dressing
37708	DraPex	45 cm x 55 cm	10 Dressing



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

TherMax



Capsicum Plaster

- Contains capsicum non-ivamide.
- New stripes coating technology allows better skin aeration.

Indications:

- Muscular pain caused by colds or strain
- Rheumatic pain
- Lumbago
- Sciatica
- Neuralgia

Directions of use

- 1) Clean & dry the painful area before applying.
- 2) Remove the protection cover of TherMax.
- 3) Apply TherMax over the affected area & press until it adheres firmly to the skin.
- 4) The duration of application should be from 2-3 days.
- 5) To remove, carefully lift up one corner of TherMax then pull the dressing off with a quick movement

Availability

Code	Description	Size	Pack
33001	Thermax Capsicum	11 cm x 28 cm	25 PCS



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

BonPlast



Elastic Adhesive Bdges

BonPlast elastic adhesive bandages are made from a high quality cotton and rayon cloth, fully spread with a stripe adhesive.

BonPlast has been designed to provide firm support and compression, for treatment of sprains and strains.

Indications:

- Varicose veins, ulcers, phlebitis and thrombosis.
- Sprains and strains
- Sports injuries
- Retaining dressings

How to use:

Apply the BonPlast bandage to the injured limb using the correct tension to suit the particular type of injury, the central yellow line is used to control the overlap and pressure.

Availability

Code	Description	Size
25001	BonPlast	5 cm x 4.5 mt
25002	BonPlast	7.5 cm x 4.5 mt
25003	BonPlast	10 cm x 4.5 mt



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

ALPHAPLAST



Adhesive Tapes

AlphaPlast Combines a well proven hypoallergenic adhesive and acetate rayon cloth (silk) to provide a light weight adhesive tape offering strength and rigidity.

AlphaPlast Is manufactured by a new technological procedures (Hot melt method) avoiding the organic solvents that are used in the conventional method leading to skin allergic reaction.

AlphaPlast Adhesive tape is available in a wide choice of widths and is suited ideally to the full range of dressing retention requirements and light strapping application.

AlphaPlast Is Tearable by hand both lengthways and widthways ensuring ease of application and time saving.


Availability

Code	Description	Size
22001	AlphaPlast	1.25 cm x 2 Yd
22002	AlphaPlast	2.5 cm x 2 Yd
22003	AlphaPlast	2.5 cm x 5 Yd
22004	AlphaPlast	5 cm x 5 Yd
22005	AlphaPlast	7.5 cm x 5 Yd
22006	AlphaPlast	10 cm x 5 Yd



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

HOSPATEK IRRIGATOR



Hospatek Irrigator

**Ideal For Personal Hygiene
Comfortable and Practical**

Hospatek irrigator Is manufactured from medical grade PVC and graduated to two liters for ease of measuring.

The Impermeable PVC is welded by the latest welding technology to ensure a leak proof seal.

Hospatek irrigator Has a tap to regulate the flow of fluid, has a hook to hang it anywhere, can be easily washed after usage and can be folded for easy storage.

Availability

Code	Description
53003	Enema Irrigator



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

BonCast



Synthetic Casting

BonCast Is a synthetic casting tape which consist of 14 needle knitted fiberglass substrate impregnated with Polyurethane resin in a shape of rolled stretches across length.

BonCast Provides smooth finish surface for comfort and safety, very conformable for easy application.

BonCast Is light in weight and strong for maximum mobility and durability.

BonCast Has high moisture vapor permeability, which reduces the risk of maceration.

BonCast Is recommended in secondary fractures and orthopedic remedial cast, and in cast bracing.

Application Procedures:

- 1- Apply under cast padding as required.
- 2- Wearing disposable gloves, open one bandage at a time and immerse in tap water 20-24 °C. Squeeze bandage whilst under water
- 3- Using even tension, wind bandage with 50% overlap.
- 4- Mould cast for 3-5 minutes (time depends on the water temperature) and smooth the surface of the cast with the flat of the hand.
- 5- Casts are load bearing in 30 minutes

Availability

Code	Description	Size	Color
90035	BonCast	5 cm x 360cm	Blue
90036	BonCast	7.5cm x 360cm	Blue
90037	BonCast	10cm x 360cm	Blue
90038	BonCast	12.5cm x 360cm	Blue



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

BonPad



BonPad Padding

BonPad Is a 100% viscose soft and absorbent natural orthopedic under cast padding.

BonPad Conforms well round body contours, and gives even coverage of the limb.

BonPad Reduces the risk of skin maceration and irritation, particularly in hot climate and patients with sensitive skin.

Indications:

For use under plaster of Paris casts for padding and protection.

How to use:

Apply over stockinet prior to the application of a casting material, use extra layers of padding to protect bony prominences.

Availability

Code	Description	Size
42012	BonPad	5 cm X 2.5 mt
42013	BonPad	7.5 cm X 2.5 mt
42014	BonPad	10 cm X 2.5 mt
42016	BonPad	15 cm X 2.5 mt
42017	BonPad	20 cm X 2.5 mt



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

GYPSONA BONFIX



P.O.P. Bdges

Plaster of Paris Bandage BP

- Fast Setting
- Minimum Plaster Loss
- Outstanding Mouldability
- Improved Handling
- Moist proof Packaging

Advantages:

Bonfix Plaster of Paris Bandage is made from a specially woven cloth, uniformly impregnated with the finest quality, fast setting plaster of Paris. The bandage is spooled on a round core for superior bandage stability. Moisture proof packaging improves the protection of the bandages.

Availability

Code	Description	Size	Pack
48001	Gypsona	7.5 cm x 2.7 mt	24 Bdg
48002	Gypsona	10 cm x 2.7 mt	24 Bdg
48003	Gypsona	15 cm x 2.7 mt	12 Bdg
48004	Gypsona	20 cm x 2.7 mt	12 Bdg
45001	Bonfix	7.5 cm x 2.7 mt	24 Bdg
45002	Bonfix	10 cm x 2.7 mt	24 Bdg
45003	Bonfix	15 cm x 2.7 mt	12 Bdg
45004	Bonfix	20 cm x 2.7 mt	12 Bdg



Afri Medical Company
 Industrial Area C3, 10th Of Ramadan City, Egypt
 Tel.: +20 554361684 - Fax: +20 222911498
 Web: www.afrimedical.com
 Web: www.facemask.afrimedical.com
 E-mail: info@afrimedical.com

URINE BAG



Hospatek Urine Bag

- Hospatek Urine Bag** Is manufactured from medical grade P. V. C which is translucent to facilitate checking and calibrated to 2000ml for ease of measuring.
- The impermeable P.V. C is welded by the latest welding technology to ensure a leakproof seal.
- Hospatek Urine Bag** Has a non-return value, which prevents the reflux of fluid from the bag and therefore reduces the danger of ascending infection.
- Does not have to be discarded when full but it is emptied when necessary by means of the drainage port at the bottom.


Availability

Code	Description
52001	2LT Urine Bag 2W - one side - singel pack
52002	2LT Urine Bag 2W - one side - 12 S



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

DRAINAGE BAG



fluid collection bag

- An Afri Medical Drainage set provides closed system for the reduction of infection risks.
- They are supplied as either EtO sterilized or as Non Sterile Bags.
- Bags are sterilizable by all known sterilization methods, including EtO, Gamma Radiation, and Steam Sterilization.



Availability

Code	Description
55042	Drainage Bag 3L - 150 micron - Without Components - Folded
55043	Drainage Bag 3L - 150 micron - With Components - Folded
55056	Drainage Bag 2LT with clamp and female DEHP Free one Hole
55049	Drainage Bag 10L - Without Components



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

BLOOD LINE



Haemodialysis Blood Line

Selection of Raw Materials complying with high Medical Standards.

Comprehensive documentation system and reserve samples to allow easy traceability. Safety is guaranteed through maintaining the residual ethylene oxide within FDA limits. Assured product safety by permanent inprocess control and 100%


leakage

A wide range of different configurations is available to suit all dialyzing machines



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

FISTULA NEEDLES



Fistula Needles

Ensured maximum blood flow with minimum resistance by a non- traumatic siliconized ultra-thin needle wall.

Better blood flow is also provided by a side hole (Back-Eye).

Easy painless penetration is guaranteed by optimum inclination of the cutting edge.

AVAILABILITY

Code	Description
74703	Fistula Needle 16G twin color coded
74704	Fistula Needle 16G twin
74710	Fistula needle 15G twin with protector
74715	Fistula Needle 17G Twin



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

VENAFLO POLYETHER SULFONE (VPES)



Haemodialysis

VENAFLO POLYETHER SULFONE DIALYZER:

- Made of an advanced synthetic polymer of excellent stability and biocompatibility.
- Membrane has excellent performance due to its high permeability for all uremia toxins.
- Dialyzers are highly stable to sterilization by all methods such AS ETO gas, gamma radiation, electron beam radiation and steam sterilization
- Dialyzers can be used to perform dialysis without the need to use systemic anticoagulants..

Venaflo Polyether Sulfone Capillary Dialyzer

Model	Surface area	Code ETO	Code Gamma	Dialyzer Housing	Dialyzer Potting Diffuser
VP 100	1.0 m ²	77715	77915	ID.30 mm	Short
VP 110	1.1 m ²	77716	77916		
VP 120	1.2 m ²	77717	77917	ID.34 mm	
VP 130	1.3 m ²	77718	77918		
VP 150	1.3 m ²	77720	77920		Long



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

FACE MASK



Disposable Face Mask

General Description:

Face masks latex free, fiberglass free, hypoallergenic, with ear loops and nose piece, Box of 50.

Intended Use:

Face mask is intended to be worn by non-medical personnel in daily interaction in public places to protect wearer from transfer of dust, fluid and particulate material.

Technical Specifications:

- Latex-free, hypoallergenic, fiberglass free and fluid resistance.
- Ultrasonic welded ear loops.
- High filtration efficiency for protection in many environments.
- Differential pressure (Pa/cm²): < 29.4.
- Low breathing resistance.
- 3-ply construction with 3 pleats of expandable folds to covers the area from the nose to the chin.
- Soft inner surface for comfort.
- Adjustable nose clip.

Donning Instructions:

- Place over nose and mouth with green side out & nose clip up.
- Secur ear loops around the ears.
- Secure nose clip.
- Secure mask under chin opening pleats.

Availability

Code	Description	Pack
36001	Face Mask	10 PCS
36003	Face Mask	5 PCS
36005	Face Mask	1 PC



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

SURGICAL MASK



Disposable Surgical Mask

General Description:

Surgical masks latex free, fiberglass free, hypoallergenic, with ear loops and nose piece, Box of 50.

Intended Use:

Surgical mask is intended to be worn by medical personnel during surgical or other medical procedures to protect both the patient and the operating personnel from transfer of microorganisms, body fluid and particulate material transfer.

Technical Specifications:

- Latex-free, hypoallergenic, fiberglass free and fluid resistance.
- Ultrasonic welded ear loops.
- High filtration efficiency for protection in many environments.
- Bacterial Filtration Efficiency (BFE): > 95%.
- Differential pressure (Pa/cm²): < 29.4.
- Low breathing resistance.
- 3-ply construction with 3 pleats of expandable folds to covers the area from the nose to the chin.
- Soft inner surface for comfort.
- Adjustable nose clip.

Donning Instructions:

- Place over nose and mouth with blue side out & nose clip up.
- Secure ear loops around the ears.
- Secure nose clip.
- Secure mask under chin opening pleats.

Availability

Code	Description	Pack
36002	Surgical Mask	10 PCS
36004	Surgical Mask	5 PCS
36006	Surgical Mask	1 PC



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

KN95 PROTECTIVE RESPIRATOR



Protective Respirator

General Description:

KN95 Respirator with Antimicrobial/Antiviral filtering facepiece respirator (FFR) for use by healthcare staff, 10 pcs per box.

Intended Use:

KN95 surgical respirator is flat folded and expands into a convex-shaped mask with elastic earloops to secure the mask to the user's face, and a strip positioned above the nose for a tighter seal around the nose and face. The device is intended for occupational use to be worn by operating room personnel during surgical procedures to help to protect both; the surgical patient and the operating room personnel from the transfer of microorganisms and body fluids. Also, to help reduce wearer exposure to pathogenic biological airborne particulates, and to protect against the transfer of particulate material.

Technical Specifications:

- Latex-free, hypoallergenic, fiberglass free and fluid resistance.
- Ultrasonic welded ear loops.
- High filtration efficiency for protection in many environments.
- Filter performance: = 95%.
- Inhalation resistance: = 343 Pa.
- Exhalation resistance: = 245 Pa.
- Multiple-ply construction with expandable fold to covers the area from the nose to the chin.
- Soft inner surface for comfort.
- Adjustable nose clip.

Donning Instructions:

- Flatten nose clip.
- Place over nose and mouth with convex side out & nose clip up.
- Secure ear loops around the ears.
- Secure nose clip.

Availability

Code	Description	Pack
36008	Kn95 Respirator	10 PCS



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com

N-SERIES FACE MASK



N-Series Facemask

General Description:

N-Series Facemask features a lightweight, 6-layer design, this face mask offers superior protection against dust and particulate material.

Intended Use:

N-Series Facemask is flat folded and expands into a convex-shaped mask with elastic ear-loops to secure the mask to the user's face, and a strip positioned above the nose for a tighter seal around the nose and face. The device is Intended to be worn by the general population for the protection from bodily fluids, droplets, and particulate material.

Technical Specifications:

- Flat fold convenience with contour molded comfort.
- Quick and easy to put on, unlike any other flat fold.
- Latex-free, hypoallergenic, fiberglass free and fluid resistance.
- Ultrasonic welded ear loops.
- Ear loops are securely attached so they don't break at the seal.
- Adjustable nose clip.
- Soft spun lining for more comfort and durability.
- Mask is naturally contoured for less pressure.
- Packed in sets of 5 for hygiene.
- Not for medical staff.

Donning Instructions:

- Flatten nose clip.
- Place over nose and mouth with convex side out & nose clip up.
- Secure ear loops around the ears.
- Secure nose clip.

Availability

Code	Description	Pack
36007	N-Series Facemask	10 PCS



Afri Medical Company
Industrial Area C3, 10th Of Ramadan City, Egypt
Tel.: +20 554361684 - Fax: +20 222911498
Web: www.afrimedical.com
Web: www.facemask.afrimedical.com
E-mail: info@afrimedical.com